

# Analysing Diffraction Data of $\text{La}_{2-x}\text{Sr}_x\text{CuO}_{4+y}$ Using McStas Virtual Experiments

Linda Udby @ Risø DTU & @ Copenhagen University


P.K. Willendrup, T.B.S. Jensen & N.H. Andersen @ Risø DTU

K. Lefmann @ Copenhagen University

Ch. Niedermayer @ PSI


B.O. Wells @ University of Connecticut

- Sr doping:  $T_M$  peaking at  $n_h \sim 1/8$ .  $T_c$  peaking at  $n_h \sim 0.16$
- (High) O doping:  $T_M = T_c = 40\text{K}$ . Staging for  $y > 0.06$
- Sr + O doping:  $T_M = T_c = 40\text{K}$ . Staging for  $x < 0.06$




H.E. Mohottala *et al*, Nature Materials **5** 377 (2006)

Undoped AFM:




LNSCO IC AFM:





N.B. Christensen et al. PRL **98** (2007)

J.M. Tranquada et al (1995-)

LSCO+O:




x	$\delta_H[\text{rlu}]$	$\delta_K[\text{rlu}]$
0	0.130(3)	0.125(3)
0.04	0.124(3)	0.124(3)
0.065	0.122(3)	0.122(3)
0.09	0.115(3)	0.123(3)





$x = 0.065$  ( $x = 0.04$ )

5% diff in proj of moment

but > 50% int diff high/low Q








$$\begin{aligned}\frac{d\sigma_M}{d\Omega}(\mathbf{Q}) &= N \frac{2\pi}{V_0} \sum_{\tau} |pf(\mathbf{Q})w^{-W}|^2 |F_{M\perp}(\mathbf{Q})|^2 \delta(\mathbf{Q} - \tau) \\ F_{M\perp}(\mathbf{Q}) &= \mathbf{Q} \times (F_M(\mathbf{Q}) \times \mathbf{Q}) \\ F_M(\mathbf{Q}) &= \sum_d \mathbf{m}_d e^{i\mathbf{Q} \cdot \mathbf{r}_d}\end{aligned}$$

- The intensity difference in high/low Q peaks is qualitatively explained by a model with stripes in 3 layers.

LCO:




LSCO+O:



• Staging number =  $1/\delta_L$

B.O. Wells *et al.*, Zeitschrift f. Physik B **100** (1996)



x	Staging number
0	3.8
0.04	5.5
0.065	disordered tilts
0.09	none

- The staging number increases with Sr content
- Decreasing O content with increasing staging number

- Intensity - cross-section
- Position - periodicity
- Width - domain size



- Intensity - cross-section
- Position - periodicity
- Width - domain size
- But often hard to know the resolution in detail in order to deconvolute domain size
- Would be nice to have resolution limited linescans for comparison

- Intensity - cross-section
- Position - periodicity
- Width - domain size
- But often hard to know the resolution in detail in order to deconvolute domain size
- Would be nice to have resolution limited linescans for comparison
- Virtual experiments on homogenous samples!

## McStas VE HOW-TO

INSTRUMENT


- Build instrument including all collimators
- Test instrument on powder 2  $\theta$  scan
- Test instrument energy resolution on Vanadium



SINGLE CRYSTAL

- Set samplesize
- Find sample mosaicity vs. rocking curve
- Find uncertainty in lattice parameter vs longitudinal scan

MORPHOLOGY

- Perform VE corresponding to diffraction data
- Deconvolute correlation length








Effective instrumental resolution and thereby the linewidth of a particular scan depends on




- Divergence of the beam before the monochromator: Size of source, geometry and m-values of guide elements
- Mosaicity of the monochromator - listed by producer as  $37'$
- Mosaicity of the analyser
- Geometrical factors: Sizes of components and distances between them
- Divergence of collimators
- Point-spread function of the position sensitive detector (PSD)
- Absolute energy of the incoming and scattered beam
- Sample parameters
  - shape/size (all samples including incoherent scatterer, powder, single-crystal)
  - particle size in sample (powder)
  - mosaicity and uncertainty in lattice parameters (single-crystal)



The samples in the virtual experiments are homogeneous and by deconvoluting the effective resolution the linewidth broadening due to finite size effect is found.

## Vanadium sample



## PowderN sample






● McStas ○ RITA-II

Fundamental reflections  
(020) and (-200)


used to set mosaicity and  $\frac{\Delta d}{d}$

Simulated (resolution) width  
agrees within 5% of measured

$x=0.09$




- The IC AFM peak is narrow and close to resolution limited...

$x=0.09$ 

- The IC AFM peak is narrow and close to resolution limited... but effects larger than 10% are observable by McStas VE



$x=0.09$



IC AFM broadened:

Fit to Voigt

$$lw=0.006(2) \text{ \AA}^{-1} \Rightarrow \text{Corr.} = 340 \pm 130 \text{ \AA}$$


Fit to Gaussian


$$w = \sqrt{w_m^2 - w_r^2}$$

$$w=0.010 \text{ \AA}^{-1} \Rightarrow \frac{2\pi}{w} \sim 610 \pm 140 \text{ \AA}$$

- The IC AFM peak is narrow and close to resolution limited... but effects larger than 10% are observable by McStas VE

$x=0.04$




$x=0.04$ 

Bmab broadened:  $w=0.034 \text{ \AA}^{-1}$ 

Domain size  $\frac{2\pi}{w} \sim 185 \pm 30 \text{ \AA}$ 

Staging broadened:  $w=0.064 \text{ \AA}^{-1}$ 

Domain size  $\frac{2\pi}{w} \sim 100 \pm 8 \text{ \AA}$

$x=0.04$ Bmab broadened:  $w=0.034 \text{ \AA}^{-1}$ Domain size  $\frac{2\pi}{w} \sim 185 \pm 30 \text{ \AA}$ Staging broadened:  $w=0.064 \text{ \AA}^{-1}$ Domain size  $\frac{2\pi}{w} \sim 100 \pm 8 \text{ \AA}$ 

- Mcstas VE is a useful tool when there is no instrumentally resolved peak in the vicinity

I would like to thank my coworkers:

- **UConn**: Barret O. Wells, Hashini E. Mohottala, Samuel B. Emery
- **McStas**: Peter K. Willendrup, Erik Knudsen, Kim Lefmann, Emmanuel Farhi
- **Risø DTU**: Niels H Andersen, Thomas B.S. Jensen, Asger B. Abrahamsen, J. i Hjøllum and Bente Lebech
- **PSI**: Niels B. Christensen, Christof Niedermayer

PLEASE VISIT OUR POSTER **I1.2** IN THE POSTER  
SESSION TUESDAY EVENING!

THANK YOU FOR LISTENING!